


# Accurate Color, Trusted Clarity for Medical Endoscopy

The MS192SA is a 19-inch endoscopy monitor tailored for digital video endoscopy and ENT applications. With a resolution of 1280×1024 and a wide color gamut, it delivers crisp, true-to-life visuals critical for clinical precision.

Equipped with Reshin's advanced color calibration and multiple GAMMA curve settings, it ensures consistent grayscale and color accuracy across diverse medical scenarios. A variety of signal interfaces allows seamless integration with a wide range of medical imaging systems.

### Accuracy, 19-inch HD images

Designed for clear visualization of 2D color video images from ENT, endoscopic, and other compatible medical imaging systems.



### **True-to-life Color Reproduction**

BT.709 Real-color Display. Enhanced by Reshin's patented technology for accurate and true-to-life color reproduction.



### **Cost-effective Option for Doctors**

Engineered with medical-grade High-bit Gamma processing and optimized gamma correction, the MS192SA ensures precise tissue differentiation while maintaining budget-friendly operation.



### **Flexible Interface Options**

It comes standard with DVI, VGA, Video, S-Video and RGB signal ports, offering seamless integration with both digital and legacy endoscopic systems. This comprehensive connectivity suite ensures compatibility without requiring additional adapters.





### **Durable Design, Rigorous Test**

Engineered for high reliability, the design has passed rigorous testing, with a substrate MTBF exceeding 50,000 hours.




## **Stand Support**

Standard Reshin customized stand with a clean, modern design. It offers height adjustability and ensures stable placement of the monitor on an endoscopy tower or table.



### **User-Friendly Control Panel and Wide Range of Connectors**

Versatile input/output options ensure seamless O.R. integration, while the user-friendly control panel allows for effortless operation.



SOURCE Hotkey to choose signal

MODE Hotkey to switch gamma mode


MENU Menu key

UP Hotkey to choose up/ Adjust contrast ratio

**DOWN** Hotkey to choose down/ Adjust brightness

**EXIT** Exit/ Auto adjustment (VGA)

POWER Power switch






# **Specifications**

19 " HD Endoscopy Monitor MS192SA		
Picture Performance		
Active Screen Size	19 inches(48.26cm)	
Pixel Pitch (H x V)	0.2928mm x 0.2928mm	
Resolution (H x V)	1280×1024	
Aspect Ratio	5:4	
Backlight	LED	
Panel Technology	TFT	
Luminance (Panel Specification)	200 cd/m² (default)	
Contrast Ratio	1000:1	
Colors Approx	16.7 billion colors	
Response Time	14ms (TYP)	
Viewing Angle (Panel Specification)	178° V. / 178° H.	
GAMMA	GAMMA1.8,GAMMA1.9,GAMMA2.0,GAMMA2.1,GAMMA2.2, GAMMA2.3,GAMMA2.4,GAMMA2.5,GAMMA2.6,DICOM,LINEAR	


### **FRONT**



#### INPUT

Digital Signal	DVI connector (x1)
Analog Signal	Video connector(x1), S-Video connector(x1), RGB/ YPbPr connector(x3), VGA connector(x1)
Sync & Control	CSYNC connector(x1), RS232 connector(x1)
USB INPUT	USB connector(x1)


## SIDE



### OUTPUT

Digital Signal	N/A
Analog Signal	$\label{lem:video} Videoconnector(x1), S-Videoconnector(x1), RGB/YPbPrconnector(\times 3)$
Sync & Control	CSYNC connector(x1)

### **BACK**



### GENERAL

02.12.012	
Power Requirements	DC 24V 5A
Power Consumption	35W (MAX)
Operating Temperature	0°C~40°C
Operating Humidity	10% to 80% (no condensation allowed)
Storage/Transport Temperature	-20°C~50°C
Storage/Transport Humidity	10% to 90% (no condensation allowed)
Dimensions (without Stand)	428.5(W)x353.5(H)x58(D)mm
Weight (without Stand)	4.8KG±0.5KG
Weight	7.3KG±0.5KG
Mounting	VESA 100x100mm
Structural Appearance	White Metal Shell
Application	ENT/ Gynecological Endoscopy/ Urology Endoscopy/ Dental Endoscopy/ Veterinary Endoscopy/ Portable Endoscopy/ Basic Laparoscopy
Supplied Accessories	LCD display (x1), Certificate (x1), Warranty card (x1), Power Cable (x1), Power Adapter (x1), Packing list (x1), DVI cable (x1), BNC cable(x1)

Shenzhen JLD Display Expert Co., Ltd all rights reserved.
Without the written permission of the copyright owner, it is not allowed to copy, distribute, delete, modify, extract, reprint or spread through networks information.
Product configuration, specifications, colors and pictures are for reference only. Subject to change without notice.
Shenzhen JLD Display Expert Co., Ltd. reserves the rights to change product design and specifications.